



Investigating instrument resolution with McStas

ISIS 2025 McStas School

PRESENTED BY MADS BERTELSEN

2025-11-05

Resolution

What do we mean by instrument resolution?

Two important concepts:

- **Point spread function:** Distribution on detector from specific $\mathbf{q}$, E
 - Used to build expected detector response from known $\mathbf{q}$, E distribution
- **Detector response:** Distribution of $\mathbf{q}$, E measured by part of detector (pixel, time)
 - Used to build up map of $\mathbf{q}$, E measured from detector response
- Depends on almost all aspects of instrument:
 - Source size/time structure, optics, sample size/shape, detector, ...



Resolution



Tools for point spread function

Spot_sample (in Contrib, by Garrett Granroth)

- Number of spots in $\mathbf{q}$, E
- Probes the point spread function for all these spots in one simulation




Resolution


Tools for detector response

- Continuous measurements **Res_sample.comp** / **Res_monitor.comp**


Start of instrument



Res_sample



Res_monitor



Sample sends rays towards monitor.
Selects random final energy in given interval.

Sample records **q**, **E** for each event
Monitor histograms those that reach it

Resolution

Tools for detector response



Terminal tool: mcresplot

```
Test_TasReso_20251031_193910 — python3 · mcresplot reso.dat — 92x61

(mcstas_fresh_2) madsbertelsen@CIN-101563 Test_TasReso_20251031_193910 % mcresplot reso.dat
This is a covariance matrix calculator using neutron events,
doi: 10.5281/zenodo.4117437,
written by T. Weber <tweber@ill.fr>, 30 March 2019.

Mean (Q, E) vector in lab system:
[-1.1213e+00  4.1128e-04  5.6444e-01 -2.7778e-03]

Covariance matrix in lab system:
[[ 9.3978e-05 -1.2078e-06  1.7440e-05  3.4412e-04]
 [-1.2078e-06  3.4789e-03  3.3121e-06  4.9043e-05]
 [ 1.7440e-05  3.3121e-06  6.1003e-05  3.3635e-04]
 [ 3.4412e-04  4.9043e-05  3.3635e-04  2.9168e-03]]
```




Resolution




Tools for detector response TOF case

- Time of flight measurements **TOFRes_sample.comp** / **TOFRes_monitor.comp**


Start of instrument



TOFRes_sample



TOFRes_monitor



Sample sends rays towards monitor.
Selects random final time in given interval.

Sample records **q**, **E** for each event
Monitor histograms those that reach it

Resolution



Tools for detector response TOF case

Terminal tool: mcresplot

```
Test_TOFRes_sample_20251031_192847 — python3 - mcresplot TOFRes.dat — 100x61
(mcstas_fresh_2) madsbertelsen@CIN-101563 Test_TOFRes_sample_20251031_192847 % mcresplot TOFRes.dat
This is a covariance matrix calculator using neutron events,
doi: 10.5281/zenodo.4117437,
written by T. Weber <tweber@ill.fr>, 30 March 2019.

Mean (Q, E) vector in lab system:
[-1.7315e+00  2.2029e-03  4.3193e+00  5.0343e+01]

Covariance matrix in lab system:
[[ 5.1688e-04  8.8460e-06  1.0220e-03  2.1371e-02]
 [ 8.8460e-06  9.7411e-04 -9.7686e-05 -1.9976e-03]
 [ 1.0220e-03 -9.7686e-05  3.5230e-03  7.3272e-02]
 [ 2.1371e-02 -1.9976e-03  7.3272e-02  4.5502e-01]]
```



Resolution

Current situation

- Separate for detector response and point spread function
 - Point spread function can be determined for n spots at once with Spot_sample


- For detector response
 - Existing solutions look at a single pixel (and time bin for TOF)
 - Getting resolution function for whole detector requires
 - Scan over instrument settings (incoming / final energy)
 - Either scan over all pixel positions or large number of components
 - For TOF instruments, scan over investigated time bin at given pixel
 - These components are admittedly designed for TAS / simple TOF instruments

Resolution

Method for full detector resolution

- Want to get both detector response and point spread function
- Want to get entire detector and all times in one simulation

Start of instrument



Before sample: $\mathbf{k}$, neutron id

After sample: $\mathbf{k}$, neutron id

Detector: pixel position, time and neutron id

Sample sends rays towards monitor.
Selects random final energy in given interval.

Resolution

Method for full detector resolution

- For each event recorded in the detector we know
 - From the detector itself: pixel position, time and neutron id
 - From this id, we can look up: initial and final wavevector -> $\mathbf{q}$, E
- With this information we can look at
 - Small subset of pixel / time positions -> Detector response for that area in $\mathbf{q}$, E
 - Small subset of $\mathbf{q}$, E -> Point spread function on detector
- This requires large dataset and code to handle data
- Still requires multiple simulations scanning over instrument settings